



# Quick Tents

4 and 9 m<sup>2</sup> Tents



Emergency

These quick assembly tents are especially designed for the care of the ill or wounded in an urban environment, ensuring privacy to the patient and the medical personnel to carry out necessary manoeuvres.

## TECHNICAL SPECIFICATIONS

Manufactured with *polyester fabric and a pre-installed fibre glass and aluminium structure.*

Deployment time is less than one minute.

Canvas: 190 g/m<sup>2</sup> 100% waterproof Polyester.

### Dimensions:

#### 4 m<sup>2</sup> Tent:

- Usable space is 1.95 x 1.95 x 1.5 m (A x P x H).
- Ø0.24 x 1.35 m dimensions transport bag (approximately).
- *Weight: 12 Kg.*

#### 9 m<sup>2</sup> Tent:

- Usable space is 3.0 x 3.0 x 1.8 m (A x P x H).
- Ø0.35 x 2.1 m dimensions transport bag (approximately).
- *Weight: 20 Kg.*