



Flexible hoses designed for the conduction of medical gases from the pipe network outlets to the points of consumption such as anaesthesia appliances, ventilators, etc...

- Available for Oxygen (white), Vacuum (yellow), Medical Air (black and white), Nitrous Oxide (blue), Carbon Dioxide (grey) and Nitrogen (black).
- Customized length and topped at the ends with the connectors specified by the client.

TECHNICAL SPECIFICATIONS

Materials:	Composed by a PVC interior pipe Middle polyester fibre reinforcement An exterior PVC cover
Antistatic:	Yes
Maximum test pressure:	15 bar
Maximum work pressure:	10 bar
Dimensions:	Vacuum: inner Ø8 mm and outer Ø14 mm Other gases: inner Ø6 mm and outer Ø12 mm
Branded with:	LOT, GAS and ISO 5359
Standard:	In compliance with UNE-EN ISO 5359