



The humidifier bottle is designed to increase the concentration of water during oxygen supply.

- Reusable after sterilization thanks to its metallic connection and the materials used in the manufacture process.
- The disposable humidifier bottle is for one use only.

TECHNICAL SPECIFICATIONS	
Humidifier bottle	
Materials:	PC, PA, Brass and NBR
Outlet connection:	DISS 9/16"
Sterilization:	Autoclave up to 121°C, 1,3 kPa (máx. 40 ciclos)
Overpressure valve:	1,4 bar
Masimum level:	200cc.
Minimum level:	100 cc.
Standards:	In compliance with UNE-EN ISO 8185
Disposable humidifier bottle	
Materiales:	Poliestireno, Poliamida, Polietileno y Silicona
Outlet connection:	DISS 9/16"
	One use only
Overpressure valve:	1,4 bar
Maximum level:	200cc
Minimum level:	100 cc
Packaged:	36 units boxes
Standards:	In compliance with UNE-EN ISO 8185

