



Accessories

Collection jars

The Collection Jar is designed to collect suction fluids. It incorporates an overflow trap and is available in different volumes: 300, 500, 1000 and 2000 ml.

- Reusable after sterilization thanks to its metallic connection and the materials used in the manufacture process.



TECHNICAL SPECIFICATIONS

300 cc jar

Materials:	Manufactured in PC, PA, chrome-plated brass and NBR
Connections:	DISS 3/4-16 UNF
Sterilization:	Sterilizable in autoclave up to 121 °C and 1,3 kPa (max. 40 cycles)
Maximum level:	11VSR0300: 300 cc.



500 cc, 1000 cc and 2000 cc jar

Materials:	PC, PA y lacquered steel (RAL 9010)
Connection:	Flexible hose
Sterilization:	Sterilizable in autoclave up to 120 °C
Maximum level:	11VSR0500: 500 cc 11VSR1000: 1000 cc 11VSR2000: 2000 cc

