

Vacuum regulators Medium Vacuum

The medium vacuum regulator **JUST VAC** provides regulation of continuous suction for its use in tracheal and pharyngeal ways during surgical procedures and continuous gastric drainage.

- The vacuum regulator provides suction from 0 to 450 mmHg.
- They are classified as **Medium Vacuum / High Flow**.
- For the safety of the patient, in case of accidental connection to a positive gas pressure, the regulators can be equipped with an optional positive pressure relief valve.
- **JUST VAC** incorporates an **ON-OFF** dispositive, through which the aspiration can be switched off without having to close the suction control valve.
- It can be provided individually or on a dual version, and it can be mounted either with a direct connector, a clamp to a rail or a Ø38.0 mm bar flange. It is supplied by a flexible hose with different European and American connectors (see *Connectors and Outlets*).
- Pressing the PUSH button, the accessories like filters, filters with overflow trap, DISS connectors or metallic connectors for flexible hoses are released for its easy replacement, ensuring the safety of the user and the patient (please see the *PUSH Accessories catalogue*).
- The original set includes the regulator, a disposable filter with safety overflow trap, three disposable filters and a metallic connector for flexible hose.



TECHNICAL SPECIFICATIONS

Materials:	ABS (yellow/black for MV), stainless steel, brass and NBR.
Modes:	ON OFF
Ranges:	Medium vacuum: 0-450 mmHg.
Supply:	from 0 to 80 L/min.
Vacuum gauge's precision:	± 2,5 % FS
Inlet connection:	1/4" BSP in stainless steel
Outlet connection:	PUSH System for a quick and easy replacement of accessories
Dimensions:	158 x 89 x 63 mm
Weight:	307 gr (without connection).
Standards:	In compliance with 93/42/CEE Directive In compliance with UNE-EN ISO 10079-3 standard
GMND Code:	44809

	DIN 13260	NF S90-116	BS 5682
Just Vac MV + Connector	11VPMVAEDN	11VPMVAENF	11VPMVAEBS
Just Vac MV Double/M + Connector	11VPMAREDN	11VPMARENF	11VPMAREBS
Just Vac MV to rail - Knob clamp	13VPMVAEPD	13VPMVAEPN	13VMAVAEPB
Just Vac MV to rail Double/M - Knob clamp	13VPMAREPD	13VPMAREPN	13VPMAREPB
Just Vac MV to Ø38,0 mm flange	14VPMVAEPD	14VPMVAEPN	14VPMVAEPB
Just Vac MV Double/M to Ø38,0 flange	14VPMAREPD	14VPMAREPN	14VPMAREPB

¹ Other outputs replace **E** (*Disposable filter*) with: **D** (*DISS output*), **T** (*filter with safety overflow trap*), **F** (*bottle of 0,5 l*), **M** (*Metallic nipple*) or **O** (*free output*).

² Other rail mounting replace **P** (*knob clamp*) with: **G** (*pushbutton clamp*), **R** (*wall plate*) or **H** (*bearing plate*).

³ Please consult our staff for other gas standars.

⁴ Please consult our staff for other models.



C/ Tungsteno 11-15
 Arganda del Rey, 28500, Madrid
 Telf.: +34 91 871 99 50
 www.itcsal.com itcsal@itcsal.com