



# Accessories

Oxygen therapy equipment



The oxygen therapy equipment is especially designed for its use in home oxygen therapy services or in hospital facilities with no oxygen piped network available.

- Gas regulator with piston flowmeter to use on an oxygen cylinder at 200 bars.
- Turning the front control, the selected flow that is supplied by the inferior connection will appear on its display.
- It incorporates a manometer to read the pressure of the cylinder.
- It also includes a humidifier bottle.

## TECHNICAL SPECIFICATIONS

- Cart of 340 x 316 x 1110 mm dimensions, manufactured in lacquered steel, with Ø200 mm wheels and security belts for the cylinder.
- Gas regulator with dial flowmeter with scales at: **0 - 0,5 - 1 - 2 - 3 - 4 - 5 - 6 - 8 - 10 - 12 and 15 l/min.** (other scales available).
- 10 l oxygen cylinder (2.000 l of oxygen).
- Humidifier bottle.
- Weight: 22,0 Kg.