

Outlet converter



Just Plug C is especially useful when the equipment connectors are different in standards from those of the valves..

- The set is composed by a connector and an outlet for the **same gas** but in compliance with **different standards**.

	SS 8752430	BS 5682	Carburos Metálicos	DIN 13260-2	Fro	Nassetti	Manufacturas Médicas	NF S90-116/FD S90-119	Oxifar
SS 8752430			✓					✓	
BS 5682								✓	
Carburos Metálicos	✓			✓	✓	✓	✓	✓	
DIN 13260-2		✓	✓	✓	✓			✓	✓
Manufacturas Médicas			✓					✓	
NF S90-116/FD S90-119	✓	✓	✓	✓	✓		✓		✓

TECHNICAL SPECIFICATIONS

Materials: Chromed-plated brass, stainless steel and NBR.

Standards: NF S 90-116, DIN 13260-2, BS 5682, SS 87 524 30, CGA V-5 (DISS)
Consult our staff about other standards