

Outlet duplicators



Just Plug D is especially designed to achieve two outlets for the same gas, in systems where there is only one.

- The set is composed by a connector, a distributor and two outlets for the same gas in compliance with the same or a different standard.

	SS 8752430	BS 5682	Carbuos Metálicos	DIN 13260-2	FTo	Manufacturas Médicas	NF S90-116/FD S90-119	Oxifar
SS 8752430	✓							
BS 5682		✓						
Carbuos Metálicos			✓				✓	✓
DIN 13260-2			✓	✓	✓		✓	
Manufacturas Médicas			✓			✓		
NF S90-116/FD S90-119	✓			✓			✓	

TECHNICAL SPECIFICATIONS

Materials: Chromed-plated brass, stainless steel and NBR.

Standards: NF S 90-116, DIN 13260-2, BS 5682, SS 87 524 30, CGA V-5 (DISS)
Consult our staff about other standards.