

Air valves for medical gases, Mini system



Specially designed sockets to have a gas outlet not anchored to the wall.

- Available for oxygen, vacuum, medicinal air and nitrous oxide (other gases consult).
- Manufactured according to European and American standards.
- Available with male thread 1/8", 1/4" (for vacuum) and with spigot for tube (please consult our staff about other gases).

Standard Type	Available	Standard Type	Available	Standard Type	Available	Standard Type	Available	Standard Type	Available	Standards Type	Available
SS 8752430	✓	ENV 737-6	AVAILABLE SOON	Nassetti	✗	NF 590-116/FD 590-119	✓	AS 2896	AVAILABLE SOON	ISO -9170-2	✗
BS 5682	✓	Fro	✗	CGA-V5 (DISS)	✓	OHMEDA	AVAILABLE SOON	UNI 9507	AVAILABLE SOON		
Carburos Metálicos	✓	Boc Saccab	AVAILABLE SOON	ISO 18082 (NIST)	AVAILABLE SOON	CHEMETRON	AVAILABLE SOON	SANS 1409	AVAILABLE SOON		
DIN 13260-2	✓	JIS T 7101	AVAILABLE SOON	Manufacturas Médicas	✓	PURITAN BENNETT	AVAILABLE SOON	Oxifar	AVAILABLE SOON		

TECHNICAL SPECIFICATIONS

Materials: Chromed-plated brass, stainless steel and NBR

Standards: NF S 90-116, DIN 13260-2, BS 5682, SS 87 524 30, CGA V-5 (DISS)
Consult our staff about other standards