

Standard connectors



Just Plug P connectors with thread for equipment or for a flexible hose.

- Available for Oxygen, Vacuum, Medical Air and Nitrous Oxide (please consult our staff about other gases).
- Manufactured according to European and American standards.
- BSP 1/4" x 19 thread for vacuum and BSP 1/8" x 28 for the other gases.

Standard Type	Available	Standard Type	Available	Standard Type	Available	Standard Type	Available	Standard Type	Available	Standards Type	Available
SS 8752430	✓	ENV 737-6	AVAILABLE SOON	Nassetti	✓	NF 590-116/FD 590-119	✓	AS 2896	AVAILABLE SOON	ISO -9170-2	✓
BS 5682	✓	Fro	✓	CGA-V5 (DISS)	✓	OHMEDA	AVAILABLE SOON	UNI 9507	AVAILABLE SOON		
Carburos Metálicos	✓	Boc Saccab	✓	ISO 18082 (NIST)	✓	CHEMETRON	AVAILABLE SOON	SANS 1409	AVAILABLE SOON		
DIN 13260-2	✓	JIS T 7101	AVAILABLE SOON	Manufacturas Médicas	✓	PURITAN BENNETT	AVAILABLE SOON	Oxifar	✓		

TECHNICAL SPECIFICATIONS

Materials: Chromed-plated brass.

Standards: NF S 90-116, DIN 13260-2, BS 5682, SS 87 524 30, CGA V-5 (DISS), ISO 18082 (NIST), UNE-EN ISO 9170-2 Type 1L SEGA