



Accessories

Infusion pumps holder

Supports specially designed for the anchoring of perfusion pumps in the patient's bedside without the need to rest on the floor, both in its rail version and its Ø38.0 mm bar version.



Rail system

TECHNICAL SPECIFICATIONS

Infusion pump holders for rail (ref.: 13SPBM0000)

Materials:	Anodized aluminium
Anchoring system:	Anchored to two rails by means of pushbutton clamps
Dimensions:	Ø25 x 1180 mm

Infusion pump holders for rail (ref.: 13SCCG3500)

Materials:	Aluminio Anodizado
Anchoring system:	Anchored to a rail by means of a clamp with knob
Dimensions:	Ø25 x 350 mm