



# Accessories

## Sphygmomanometers to rail

Aneroid mercury sphygmomanometers, with sleeve for adults or children, for hospital use.

- Adaptable to rail by clamps. This enables its mobility between patients.
- It can also be mounted on a wall.



### TECHNICAL SPECIFICATIONS

- Round body of resistant and sturdy plastic, 65° steerable to each side.
- Copper-beryllium membrane.
- No zero point fixation required in  $\varnothing 147$  mm with scale up to 300 mmHg.
- Coil tubing extendible up to 3 m, latex bulb and chrome-plated valve.