



BASIC modular

Basic horizontal bedhead



The **Modular BASIC Bedhead unit** is a space-designed system for Hospitalization Areas. Its modular design allows several configurations that adapt perfectly to the functional and aesthetic needs of any areas of the Hospital.

Design with option of wall or wall fixing, with one or more profiles, which allows any standard or special configuration.

Easy to clean due to its straight lines.

The design guarantees an internal separation and for individualized maintenance of all the systems that make up these bedheads, both gas and electric, and has expanded its capacity to include direct and indirect lights and rail

The installation has been simplified to be carried out by a single person (always depending on the size of the headboard).

TECHNICAL SPECIFICATIONS

Materials:	Structures and covers made of extruded aluminum Covers end in ABS
Dimensions:	Without lights: 64 mm deep by 106 mm by compartments With lights: 98 mm deep by 282 mm with a compartment, increasing 106 mm per additional compartment
Finished:	Laquered main body in RAL 9010 Laquered covers in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
Standards:	In compliance with UNE-EN ISO 11197