



TECHNICAL WAY

Basic horizontal bedhead

The **TECHNICAL WAY Bedhead unit** is a wall bedhead that integrates the supply of medical gases, power, telephone and computer plugs, etc...

It can be mounted either in single, double or triple bars, effectively separating the electricity cables and the gas pipes. Our technical office will adjust the dimensions, colours and final equipment to the specific needs of each project upon demand.



Bedhead units

TECHNICAL SPECIFICATIONS

Materials:	Structure and covers manufactured in extruded aluminium Cover ends manufactured in polystyrene
Dimensions:	The dimensions are 74 mm depth by 120 mm canal
Treatments:	Main body lacquered: RAL 9010 Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
Standards:	In compliance with UNE-EN ISO 11197