



# BASIC WALL

Decorative bedhead



The BASIC Wall bedhead unit is a hospital bedhead designed to respond to the highest aesthetics and functionality requirements, becoming another architectural element of the environment.

Can be easily combined with wall light units to offer a discreet and elegant solution for normal care rooms.

The sliding frame mounted on a ball bearing system makes it easy to open for access to de operation devices, which are protected from dust.

Can be provided with indirect LED light.

**"Basic Wall is fully customized with the decoration of your choice"**

Bedhead units

## TECHNICAL SPECIFICATIONS

<b>Materials:</b>	Structure and covers manufactured in extruded aluminium and lacquered steel Cover manufactured in wood
<b>Dimensions:</b>	950 x 625 x 140 mm, protrudes from the wall surface 52 mm
<b>Finishings:</b>	Main body lacquered: RAL 9010. Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035 The front cover can be customized with different bedhead finishings for different services areas, and the image of the hospital can be incorporated
<b>Standards:</b>	In compliance with UNE-EN ISO 11197



C/ Tungsteno 11-15  
Arganda del Rey, 28500, Madrid  
Telf.: +34 91 871 99 50  
www.itcsal.com itcsal@itcsal.com

