



CLAIRE

Decorative horizontal bedhead

The **CLAIRE bedhead unit** is a hospital bedhead designed to respond to the highest aesthetics and functionality requirements, becoming another architectural element of the environment. The mechanisms, except the gas outlets, are incorporated to its lower profile, parallel to the ground, thus being hidden by the front cover.

Starting on some basic components, our technical office will adapt the dimensions, colours and final equipment to the specific needs of each project upon demand.

"New version with less depth"



Bedhead units

TECHNICAL SPECIFICATIONS	
Materials:	Structure and covers manufactured in extruded aluminium Translucent covers in Polycarbonate Covers ends in ABS Its new dimensions are 102 mm depth and 233 mm height
Finishings:	Main body lacquered: RAL 9010 Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035 The front cover can be customized with different bedhead finishings for different services areas, and the image of the hospital can be incorporated.
Standards:	In compliance with UNE-EN ISO 11197