



CLASSIC

Horizontal bedhead with lighting

CLASSIC is a wall bedhead able to integrate the supply of medical gases, electrical, telephone and computer plugs and the call device for nursing assistance, as well as to cast indirect light to the patient's bed.

With basic components, our technical office will adapt the dimensions, colours and final equipment to the specific needs of each project upon demand.



TECHNICAL SPECIFICATIONS

Materials:	Structure and covers manufactured in extruded aluminium Translucent covers in polycarbonate Cover ends in ABS
Dimensions:	117 mm depth by 238 mm height
Finishing:	Main body lacquered: RAL 9010. Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
Standards:	In compliance with UNE-EN ISO 11197