



# Accessories

## NEXUS cart - Metallic drawer

Set of drawer furniture that is attached to different structures of carts and bedheads.

- It is composed of an upper shelf and drawers in different formats and side flages that allow their installation on the carriage structure and subsequent regulation.
- The inner buckets of the drawers can easily be completely extracted out of the furniture for their cleaning or manipulation.
- Thanks to the slide partitions, it is possible to establish different divisions in the interior of the drawers.



### TECHNICAL SPECIFICATIONS

**Materials:**

Upper shelf manufactured with a Lacquered Steel Sheet  
 Drawers manufactured in ABS and available in 400x300 mm and 600x400 mm with 50, 100 and 200 in height  
 The interior dividers are in ABS  
 The exterior structure is of Extruded Aluminium Lacquered in white RAL 9010  
 Finishing elements in ABS with different colours.

**Finishing:**

Lacquered: RAL 9010  
 Finishing elements in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035