



Accessories

NEXUS twistable cart / twistable and sliding cart

This cart is capable of turning gently on a vertical axis and moving on the horizontal axis, holding equipment and accessories on the patient's bedhead, without the need of supports on the floor or on the walls.

- There are two standard lengths adaptable to each need. Both racks have four lateral profiles and two inferior ones for the anchoring of various accessories from our rail commercial line.
- Accessories such as shelves, drawers and auxiliary rails can be easily accommodated to vertical tubes with removable flanges, allowing a quick change in the configuration of the cart.



Bedhead units

TECHNICAL SPECIFICATIONS

Materials:	Lacquered steel profiles and $\varnothing 38$ mm stainless steel tubes They have a frontal and an inferior system that enables displacement and turn blocking. Available in two widths and two standard heights, although the latter can be adapted to the client's needs.
Finishing:	Lacquered: RAL 9010 Finishing elements in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
Maximum load:	150 kg