



# NEXUS LM and LX

Pendant bedhead

**NEXUS II LM and LX pendant bedheads** are designed to integrate the supply of medical gases and electricity and communications plugs, and are specially indicated for emergency wards, reanimation areas and intensive care units.

- **NEXUS II LM bedhead** has two brackets attached to the ceiling (three for the **LX** version) and a horizontal profile containing four isolated compartments, which can be equipped differently depending on the needs of each project.
- Complemented with only rotating or rotating and sliding racks of different heights, trays, drawers and diverse accessories from our rail commercial line.



## TECHNICAL SPECIFICATIONS

<b>Materials:</b>	Structure and covers manufactured in extruded aluminium Translucent covers in polycarbonate Cover ends in polystyrene
<b>Dimensions:</b>	503 mm depth by 181 mm height and lengths of 3m for the LM version and 4,5m for the LX version, being able to install units of higher length by joining several bedheads
<b>Finishing:</b>	Main body lacquered RAL 9010. Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
<b>Standards:</b>	In compliance with UNE-EN ISO 11197