



# NEXUS LT

Ceiling-pendant

The **ceiling-pendant NEXUS LT** are designed to integrate the supply of medical gases and electricity and communications plugs, and are specially indicated for surgical areas and intensive care and neonatal units.

- The ceiling-pendant comprises a support to the ceiling and a horizontal profile containing six isolated compartments (three in the LT-S version), which can be equipped differently depending on the needs of each project.
- Its design displays the gas outlets in working arrangement, allowing the operator to connect flowmeters and vacuum regulators directly during delicate situations, avoiding connection by flexible pipes.
- Complemented with rotating racks of different heights, trays, drawers and diverse accessories from our rail commercial line.



Bedhead units

## TECHNICAL SPECIFICATIONS

<b>Materials:</b>	Structure and covers in extruded aluminium Translucent covers in polycarbonate Covers ends in polystyrene
<b>Dimensions:</b>	414 mm depth by 255 mm height and 500 mm length
<b>Finishings:</b>	Main body lacquered: RAL 9010 Covers lacquered in: Red RAL 3003, Blue RAL 5015, Green RAL 6019, Grey RAL 7035
<b>Standards:</b>	In compliance with UNE-EN ISO 11197