



Stretchers and stands

Stretcher stands



Bellow' stand. Tool-free opening and closing. Blockage performed by the own stretcher.


TECHNICAL SPECIFICATIONS

Manufactured with *anodized aluminium frames*.

Dimensions:

- *Closed:* 60 x 19 x 88 cm
- *Open:* 60 x 114 x 88 cm
- *Weight:* 6,8 Kg

Pre-assembled stand. Set composed by two half stands that must be assembled before their.



TECHNICAL SPECIFICATIONS

Manufactured with *anodized aluminium frames*.

Dimensions:

- *Closed:* 6 x 59 x 144 cm
- *Open:* 120 x 59 x 78 cm
- *Weight:* 6,25 Kg



Foldable stand. A pin allows the stretcher to be folded until all components are placed on the same level.

TECHNICAL SPECIFICATIONS

Manufactured with *anodized aluminium frame*.

Dimensions:

- *Closed:* 6 x 57 x 183 cm
- *Open:* 102 x 57 x 64 cm
- *Weight:* 9,0 Kg