



PMA accessories

Oxygen equipment for field Hospitals



This distribution system allows the user to deploy a flexible oxygen tube net inside a field hospital, in order to supply specific oxygen flows to each patient.

- The rack cylinder trolley is manufactured in aluminium and allows the storage of four cylinders.
- Systems available for 4, 6 and 8 patients.

TECHNICAL SPECIFICATIONS

The Oxygen equipment for field hospitals is composed of:

- Trolley for the transport of 4 cylinders (Not included).
- Regulator equipment for 2 cylinders.
- Oxygen distribution network with polyester bridles to hang valves and 4, 6 or 8 dial flowmeters with 0 to 15 l/min flow.
- Trolley case in black plastic with wheels and a handle for the transport of all the necessary material.