



# REANIVAC

Resuscitation case



Emergency

The **REANIVAC** emergency equipments have been especially designed to be extraordinarily flexible, autonomous and easy-to-use.

**REANIVAC I** allows:

- Oxygen inhalation for recovering patients or to transfer reanimated patients.
- Oxygen enrichment of the manual resuscitator.

**REANIVAC II** also allows:

- Venturi-effect inhalations through the use of an oxygen bottle as negative pressure source.

## TECHNICAL SPECIFICATIONS

**REANIVAC I** is composed by:

- Double fold polyethylene case with 540 x 410 x 155 mm dimensions. It has an interior pocket with a 2,25 l capacity, accessible from two transparent polycarbonate easy open lids.
- **UNIVEN** with overpressure limiting valve, adult mask, extension lead for connection to the flowmeter and reservoir bag.
- 400 l oxygen cylinder with a regulator and a digital flowmeter from 0 to 15 l/min.
- Medium concentration mask for adults and children.
- Case with mouth to mouth mask.
- Set of 70, 80 and 110 mm Guedel airways.
- Instructions for use.

**REANIVAC II** is composed of every **REANIVAC I**'s material plus:

- Venturi-effect suction system, with collection jar and suction hose.
- Gold or silver thermal blanket.
- Tongue clamp.
- Jaw-spreader.