



# Accessories

On demand valve



Emergency

The on-demand valve has been designed for the patient-controlled administration of 100% oxygen.

It is specially designed for emergency services. In addition, it has a manual ventilation button located on the front of the device, which helps keep the patient's airway open.

The on-demand valve is designed to oppose very little resistance to flow, which in turn means less effort for the patient to breathe.

## TECHNICAL SPECIFICATIONS

- **Manual flow:** 40,0 l/min.
- **Flow under inhalation demand:** 0 - 160 l/min.
- **Pressure releasing valve:** 60 cm H<sub>2</sub>O with 60 l/min.
- **Resistance to exhalation:** -6,0 cm H<sub>2</sub>O with 60 l/min.
- **Inlet connection:** Quick plug NODAIR.
- Yellow spiral PA tube Ø6,0 x Ø4,0 mm.
- **Dimensions:** 9,0 x 5,0 cm
- **Weight:** 280 grs.