



Accessories

PEEP Valve



PEEP valves are designed to complement manual resuscitators, facilitating continuous pressure regulation.

- PEEP valve for adults can be set between 5 and 20 cm of H₂O.
- PEEP valve for children can be set between 1 and 10 of H₂O.

According to the Guidelines and Recommendations for Resuscitation 2015, the use of a PEEP to facilitate continuous pressure regulation is recommended in adults and especially in neonates:

"All term or preterm infants who remain in apnea despite initial measures should receive positive pressure ventilation after initial lung insufflation. A PEEP of approximately 5 cm H₂O pressure will be applied to preterm infants receiving positive pressure ventilation."

(Wyllie J, Perlman JM, Kattwinkel J, et al. Part 7: Neonatal resuscitation: 2015)

TECHNICAL SPECIFICATIONS

- PEEP valve for adults: 5-20 cm H₂O pressure.
- PEEP valve for children: 0-10 cm H₂O pressure.
 - Material: Polysulfone.
 - Connector: 18/22 mm.
 - Sterilization: Autoclave