



Stretchers and stands

Double fold stretcher



The **Double fold** stretcher is a light and rigid platform used to carry the wounded and the ill from inaccessible places to field hospitals or emergency vehicles, in order to be taken care of in a suitable environment.

- Given its short dimensions once folded in its cover, it is specially indicated both for places with short room available, and for public places where it may be seldom used.

TECHNICAL SPECIFICATIONS

Composed by a metallic articulated enclosure of 2,02 x 0,54 m dimension with four lower legs of 13 cm, handles on the edges and a canvas attached to the enclosure by means of two welded riveted openings.

The different parts that compose the structure are manufactured with light alloy and stainless steel, and are treated against corrosion.

The canvas (100 % PES PVC 900 g/m², Fireproof Cl.2) provides great resistance.

The transport bag is made of 87% PVC y 13% Poliéster.

The stackable stretcher must be cleaned with non abrasive soaps and plenty of water after its use, and then allowed to dry in a well ventilated area before its storage.

Total weight: 6,2 Kg